


LUER BONE RONGEUR STRAIGHT 17 CM 8 MM

Réf. IP48017



LUER BONE RONGEUR STRAIGHT 17 CM 8 MM



Réf. IP48017

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040026992
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,0	Width packaging (cm):	
Length (cm) :	17,0	Length packaging (cm) :	
Height (cm) :	1,4	Height packaging (cm) :	
Weight (kg) :	0,147	Weight packaging (kg) :	0,147
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:21 - This article is available on our website : www.gsh-med.fr with the following reference IP48017

All characteristics are given for information only and may be changed