

BEYER BONE RONGEUR DOUBLE ARTICULATION LIGHTLY CURVED 18 CM 3.5 MM

Réf. IP48018



BEYER BONE RONGEUR DOUBLE ARTICULATION LIGHTLY CURVED 18 CM 3.5 MM



Réf. **IP48018**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040027005
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,0	Width packaging (cm):	
Length (cm) :	18,7	Length packaging (cm) :	
Height (cm) :	1,2	Height packaging (cm) :	
Weight (kg) :	0,146	Weight packaging (kg) :	0,146
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/06/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference IP48018

All characteristics are given for information only and may be changed