

DOYEN INTESTINAL CLAMP 23 CM STRAIGHT ATRAUMATIC JAWS

Réf. IP14428



DOYEN INTESTINAL CLAMP 23 CM STRAIGHT ATRAUMATIC JAWS



Réf. IP14428

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040086019
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,2	Width packaging (cm):	
Length (cm) :	24,0	Length packaging (cm) :	
Height (cm) :	0,2	Height packaging (cm) :	
Weight (kg) :	0,066	Weight packaging (kg) :	0,066
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 17:02 - This article is available on our website : www.gsh-med.fr with the following reference IP14428

All characteristics are given for information only and may be changed