

COLLIN MALLET DOUBLE HEADS 20 CM 500 GR (STAINLESS STEEL)

Réf. IM01250



COLLIN MALLET DOUBLE HEADS 20 CM 500 GR (STAINLESS STEEL)



Réf. IM01250

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040020389
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,7	Width packaging (cm):	6,7
Length (cm) :	21,5	Length packaging (cm) :	21,5
Height (cm) :	4,0	Height packaging (cm) :	4,0
Weight (kg) :	0,637	Weight packaging (kg) :	0,637
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:49 - This article is available on our website : www.gsh-med.fr with the following reference IM01250

All characteristics are given for information only and may be changed