

DOYEN INTESTINAL CLAMP ATRAUMATIC JAWS STRAIGHT 13 CM

Réf. IP14430



DOYEN INTESTINAL CLAMP ATRAUMATIC JAWS STRAIGHT 13 CM



Réf. **IP14430**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040086026
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,5	Width packaging (cm):	
Length (cm) :	21,0	Length packaging (cm) :	
Height (cm) :	0,3	Height packaging (cm) :	
Weight (kg) :	0,025	Weight packaging (kg) :	0,025
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 17:02 - This article is available on our website : www.gsh-med.fr with the following reference IP14430

All characteristics are given for information only and may be changed