

HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°1/2

Réf. IS01401

Exemple de bougies Hégar double, inox.
15 tailles disponibles, nous consulter.



HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°1/2

Réf. IS01401

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029542
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,3	Width packaging (cm):	
Length (cm) :	19,0	Length packaging (cm) :	
Height (cm) :	0,2	Height packaging (cm) :	
Weight (kg) :	0,007	Weight packaging (kg) :	0,007
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:34 - This article is available on our website : www.gsh-med.fr with the following reference IS01401

All characteristics are given for information only and may be changed