

HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°13/14

Réf. IS01413



HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°13/14



Réf. IS01413

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029603
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,2	Width packaging (cm):	
Length (cm) :	20,3	Length packaging (cm) :	
Height (cm) :	1,2	Height packaging (cm) :	
Weight (kg) :	0,063	Weight packaging (kg) :	0,063
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:35 - This article is available on our website : www.gsh-med.fr with the following reference IS01413

All characteristics are given for information only and may be changed