

## HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°15/16

Réf. IS01415

Exemple de bougies Hégar double, inox.  
15 tailles disponibles, nous consulter.



# HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°15/16



Réf. IS01415

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029610
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,5	Width packaging (cm):	
Length (cm) :	21,0	Length packaging (cm) :	
Height (cm) :	5,0	Height packaging (cm) :	
Weight (kg) :	0,077	Weight packaging (kg) :	0,077
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:35 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IS01415

All characteristics are given for information only and may be changed