

## HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°3

Réf. IS01603



Exemple de bougies Hégar simple, inox.  
30 tailles disponibles, nous consulter



# HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°3



Réf. IS01603

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029733
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,3	Width packaging (cm):	
Length (cm) :	19,1	Length packaging (cm) :	
Height (cm) :	0,3	Height packaging (cm) :	
Weight (kg) :	0,012	Weight packaging (kg) :	0,012
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:36 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IS01603

All characteristics are given for information only and may be changed