

DESMARRES LID RETRACTOR 8 MM

Réf. IR10008



DESMARRES LID RETRACTOR 8 MM



Réf. IR10008

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040150529
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	1,0
Length (cm) :		Length packaging (cm) :	13,5
Height (cm) :		Height packaging (cm) :	1,2
Weight (kg) :	0,026	Weight packaging (kg) :	0,026
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:22 - This article is available on our website : www.gsh-med.fr with the following reference IR10008

All characteristics are given for information only and may be changed