

SAUTTER HYDRODISSECTION CANNULA

Réf. IC92118



SAUTTER HYDRODISSECTION CANNULA



Réf. IC92118

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | No | | |
| Marquage | CE | Barcode | 3661040145006 |
| Guarantee (month): | 24 | Customs Code | 901850 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 6,7 | Width packaging (cm): | 6,7 |
| Length (cm) : | 20,2 | Length packaging (cm) : | 20,2 |
| Height (cm) : | 1,0 | Height packaging (cm) : | 1,0 |
| Weight (kg) : | 0,200 | Weight packaging (kg) : | 0,014 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



19/05/2022 21:22 - This article is available on our website : www.gsh-med.fr with the following reference IC92118

All characteristics are given for information only and may be changed