

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°17

Réf. IS01617



Exemple de bougies Hégar simple, inox.
30 tailles disponibles, nous consulter



HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°17



Réf. IS01617

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe Ir | |
| Marquage | CE | Barcode | 3661040029870 |
| Guarantee (month): | 24 | Customs Code | 9018.9090 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 1,6 | Width packaging (cm): | |
| Length (cm) : | 19,1 | Length packaging (cm) : | |
| Height (cm) : | 1,6 | Height packaging (cm) : | |
| Weight (kg) : | 0,057 | Weight packaging (kg) : | 0,057 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 21:36 - This article is available on our website : www.gsh-med.fr with the following reference IS01617

All characteristics are given for information only and may be changed