

## HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°20

Réf. IS01620



Exemple de bougies Hégar simple, inox.  
30 tailles disponibles, nous consulter



# HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°20



Réf. IS01620

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029900
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,1	Width packaging (cm):	
Length (cm) :	19,0	Length packaging (cm) :	
Height (cm) :	1,8	Height packaging (cm) :	
Weight (kg) :	0,122	Weight packaging (kg) :	0,122
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:37 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IS01620

All characteristics are given for information only and may be changed