

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°21

Réf. IS01621



Exemple de bougies Hégar simple, inox.
30 tailles disponibles, nous consulter



HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°21



Réf. IS01621

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029917
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,3	Width packaging (cm):	
Length (cm) :	19,0	Length packaging (cm) :	
Height (cm) :	1,5	Height packaging (cm) :	
Weight (kg) :	0,109	Weight packaging (kg) :	0,110
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:37 - This article is available on our website : www.gsh-med.fr with the following reference IS01621

All characteristics are given for information only and may be changed