

## PROBE WITH ONE KNOB 16 CM

Réf. IS08416



## PROBE WITH ONE KNOB 16 CM



Réf. IS08416

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040118543
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,003	Weight packaging (kg) :	0,003
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:42 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IS08416

All characteristics are given for information only and may be changed