

DUPUY-DUTEMPS RASPATORY

Réf. IR09300



DUPUY-DUTEMPS RASPARTORY



Réf. IR09300

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040029399
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,9	Width packaging (cm):	
Length (cm) :	14,2	Length packaging (cm) :	
Height (cm) :	0,8	Height packaging (cm) :	
Weight (kg) :	0,022	Weight packaging (kg) :	0,022
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference IR09300

All characteristics are given for information only and may be changed