

## ROCHESTER-PEAN FORCEPS STRAIGHT 14 CM FOR LEFT HAND

Réf. IP02214G



# ROCHESTER-PEAN FORCEPS STRAIGHT 14 CM FOR LEFT HAND



Réf. IP02214G

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040041056
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :	14,0	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:57 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP02214G

All characteristics are given for information only and may be changed