

NAIL CUTTER 13 CM (SECATEUR)

Réf. IP05423



NAIL CUTTER 13 CM (SECATEUR)



Réf. IP05423

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.
 The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.
 Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.
 100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

FEATURES			
Medical device	Yes	(Sans)	
Marquage	AUCUN	Barcode	3661040021959
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,5	Width packaging (cm):	18,0
Length (cm) :	13,0	Length packaging (cm) :	6,0
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,102	Weight packaging (kg) :	0,140
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:10 - This article is available on our website : www.gsh-med.fr with the following reference IP05423

All characteristics are given for information only and may be changed