

KOCHER FORCEPS STRAIGHT 1X2 TH 14 CM (COLLIN LOCK)

Réf. IP02414



KOCHER FORCEPS STRAIGHT 1X2 TH 14 CM (COLLIN LOCK)



Réf. **IP02414**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040021690
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,0	Width packaging (cm):	
Length (cm) :	14,0	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,037	Weight packaging (kg) :	0,037
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:58 - This article is available on our website : www.gsh-med.fr with the following reference IP02414

All characteristics are given for information only and may be changed