

KOCHER FORCEPS STRAIGHT 1X2 TH 18 CM (COLLIN LOCK)

Réf. IP02418



KOCHER FORCEPS STRAIGHT 1X2 TH 18 CM (COLLIN LOCK)



Réf. **IP02418**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040106175
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,037	Weight packaging (kg) :	0,037
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:58 - This article is available on our website : www.gsh-med.fr with the following reference IP02418

All characteristics are given for information only and may be changed