

LONDON COLLEGE TWEEZERS 15 CM

Réf. IP04715



LONDON COLLEGE TWEEZERS 15 CM



Réf. **IP04715**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040021874
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,0	Width packaging (cm):	
Length (cm) :	15,0	Length packaging (cm) :	
Height (cm) :	0,6	Height packaging (cm) :	
Weight (kg) :	0,023	Weight packaging (kg) :	0,023
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:59 - This article is available on our website : www.gsh-med.fr with the following reference IP04715

All characteristics are given for information only and may be changed