

WEIGHT 400 G FOR SCHERBACK HANDLE/VAGINAL SPECULUM

Réf. IP05400



WEIGHT 400 G FOR SCHERBACK HANDLE/VAGINAL SPECULUM

Réf. IP05400

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040121321
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	9,0
Length (cm) :		Length packaging (cm) :	18,5
Height (cm) :		Height packaging (cm) :	4,0
Weight (kg) :	0,077	Weight packaging (kg) :	0,077
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 17:00 - This article is available on our website : www.gsh-med.fr with the following reference IP05400

All characteristics are given for information only and may be changed