

## MICRO HALSTEAD FORCEPS STRAIGHT 1X2 TEETH 10 CM

Réf. IP08010



# MICRO HALSTEAD FORCEPS STRAIGHT 1X2 TEETH 10 CM



Réf. **IP08010**

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.  
The steels come into contact with human tissues and must therefore comply with international standards and regulations in force. Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040022420
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,5	Width packaging (cm):	
Length (cm) :	10,0	Length packaging (cm) :	
Height (cm) :	0,3	Height packaging (cm) :	
Weight (kg) :	0,019	Weight packaging (kg) :	0,019
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 17:00 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP08010

All characteristics are given for information only and may be changed