

DOYEN NEEDLE HOLDER 13 CM

Réf. IP30413



DOYEN NEEDLE HOLDER 13 CM



Réf. **IP30413**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040024813
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :	13,0	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,038	Weight packaging (kg) :	0,038
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/11/2022 12:43 - This article is available on our website : www.gsh-med.fr with the following reference IP30413

All characteristics are given for information only and may be changed