

DOYEN NEEDLE HOLDER 15 CM

Réf. IP30415



DOYEN NEEDLE HOLDER 15 CM



Réf. IP30415

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040024837
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :	15,0	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,039	Weight packaging (kg) :	0,039
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/11/2022 12:43 - This article is available on our website : www.gsh-med.fr with the following reference IP30415

All characteristics are given for information only and may be changed