

TOHOKU TOWEL FORCEPS 10 CM

Réf. IP10615



 **Holtex**



TOHOKU TOWEL FORCEPS 10 CM



Réf. **IP10615**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040085890
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,0	Width packaging (cm):	
Length (cm) :	20,5	Length packaging (cm) :	
Height (cm) :	1,3	Height packaging (cm) :	
Weight (kg) :	0,023	Weight packaging (kg) :	0,023
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 17:01 - This article is available on our website : www.gsh-med.fr with the following reference IP10615

All characteristics are given for information only and may be changed