

LANGE-HOHMANN BONE LEVER 27 CM - 43 MM

Réf. IL11843



LANGE-HOHMANN BONE LEVER 27 CM - 43 MM



Réf. **IL11843**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040020303
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,4	Width packaging (cm):	4,4
Length (cm) :	26,7	Length packaging (cm) :	26,7
Height (cm) :	0,4	Height packaging (cm) :	0,4
Weight (kg) :	0,135	Weight packaging (kg) :	0,135
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:48 - This article is available on our website : www.gsh-med.fr with the following reference IL11843

All characteristics are given for information only and may be changed