

AIR & LIQUID INJECTION CANNULA ANGULATED 0.5 MM

Réf. IC92105



AIR & LIQUID INJECTION CANNULA ANGULATED 0.5 MM



Réf. **IC92105**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040084107
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,4	Width packaging (cm):	
Length (cm) :	5,8	Length packaging (cm) :	
Height (cm) :	1,4	Height packaging (cm) :	
Weight (kg) :	0,015	Weight packaging (kg) :	0,015
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:58 - This article is available on our website : www.gsh-med.fr with the following reference IC92105

All characteristics are given for information only and may be changed