

BENNETT BONE LEVER 24 CM

Réf. IL11924



BENNETT BONE LEVER 24 CM



Réf. IL11924

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040020327
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,2	Width packaging (cm):	4,2
Length (cm) :	24,2	Length packaging (cm) :	24,2
Height (cm) :	6,4	Height packaging (cm) :	6,4
Weight (kg) :	0,171	Weight packaging (kg) :	0,171
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/06/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference IL11924

All characteristics are given for information only and may be changed