

CORWIN TONSIL FORCEPS

Réf. IC99300



CORWIN TONSIL FORCEPS



Réf. IC99300

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040013244
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,7	Width packaging (cm):	
Length (cm) :	16,6	Length packaging (cm) :	
Height (cm) :	2,5	Height packaging (cm) :	
Weight (kg) :	0,083	Weight packaging (kg) :	0,083
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



23/11/2022 21:09 - This article is available on our website : www.gsh-med.fr with the following reference IC99300