

COLLIN MALLET DOUBLE HEADS 20 CM 200 GR (STAINLESS STEEL)

Réf. IM01220



COLLIN MALLET DOUBLE HEADS 20 CM 200 GR (STAINLESS STEEL)



Réf. IM01220

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040020372
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,4	Width packaging (cm):	4,9
Length (cm) :	20,7	Length packaging (cm) :	20,7
Height (cm) :	2,8	Height packaging (cm) :	2,8
Weight (kg) :	0,026	Weight packaging (kg) :	0,026
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:49 - This article is available on our website : www.gsh-med.fr with the following reference IM01220

All characteristics are given for information only and may be changed