

DOYEN Mallet SINGLE HEAD 21 CM 200 GR (STAINLESS STEEL)

Réf. IM02320



DOYEN Mallet SINGLE HEAD 21 CM 200 GR (STAINLESS STEEL)



Réf. **IM02320**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040020396
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,5	Width packaging (cm):	4,5
Length (cm) :	21,0	Length packaging (cm) :	21,0
Height (cm) :	3,1	Height packaging (cm) :	3,1
Weight (kg) :	0,280	Weight packaging (kg) :	0,280
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:49 - This article is available on our website : www.gsh-med.fr with the following reference IM02320

All characteristics are given for information only and may be changed