


DUPUY-DUTEMPS LACHRYMAL SAC & DACRYO RETRACTOR

Réf. IE41180



DUPUY-DUTEMPS LACHRYMAL SAC & DACRYO RETRACTOR



Réf. IE41180

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040017778
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,0	Width packaging (cm):	
Length (cm) :	7,5	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,005	Weight packaging (kg) :	0,005
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:59 - This article is available on our website : www.gsh-med.fr with the following reference IE41180

All characteristics are given for information only and may be changed