

## FINOCHIETTO NEDDLE HOLDER STRAIGHT 15 CM

Réf. IP39415



# FINOCHIETTO NEDDLE HOLDER STRAIGHT 15 CM



Réf. IP39415

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040025285
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	30,0	Width packaging (cm):	
Length (cm) :	15,2	Length packaging (cm) :	
Height (cm) :	0,6	Height packaging (cm) :	
Weight (kg) :	0,045	Weight packaging (kg) :	0,045
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 17:20 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP39415

All characteristics are given for information only and may be changed