

DISSECTOR DOUBLE END 18 CM - 4 + 5 MM

Réf. ID29045



DISSECTOR DOUBLE END 18 CM - 4 + 5 MM



Réf. ID29045

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040084374
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,037	Weight packaging (kg) :	0,037
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:29 - This article is available on our website : www.gsh-med.fr with the following reference ID29045

All characteristics are given for information only and may be changed