

ADSON ELEVATOR 17 CM

Réf. IE32217



ADSON ELEVATOR 17 CM



Réf. IE32217

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040017310
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,7	Width packaging (cm):	
Length (cm) :	16,8	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,039	Weight packaging (kg) :	0,039
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:39 - This article is available on our website : www.gsh-med.fr with the following reference IE32217

All characteristics are given for information only and may be changed