

TRACHEAL BLUNT RETRACTOR 2 TEETH

Réf. IE33216



TRACHEAL BLUNT RETRACTOR 2 TEETH



Réf. IE33216

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040017389
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,7	Width packaging (cm):	
Length (cm) :	16,0	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,022	Weight packaging (kg) :	0,022
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:39 - This article is available on our website : www.gsh-med.fr with the following reference IE33216

All characteristics are given for information only and may be changed