

## TRACHEAL BLUNT RETRACTOR 3 TEETH

Réf. IE33316



## TRACHEAL BLUNT RETRACTOR 3 TEETH



Réf. IE33316

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040017396
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,0	Width packaging (cm):	
Length (cm) :	18,2	Length packaging (cm) :	
Height (cm) :	0,3	Height packaging (cm) :	
Weight (kg) :	0,025	Weight packaging (kg) :	0,025
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



11/08/2022 21:07 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IE33316

All characteristics are given for information only and may be changed