

BAKES DUCT DILATOR 8 MM

Réf. IE39108



BAKES DUCT DILATOR 8 MM



Réf. IE39108

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040017600
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,0	Width packaging (cm):	
Length (cm) :	30,9	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,032	Weight packaging (kg) :	0,032
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



23/11/2022 21:10 - This article is available on our website : www.gsh-med.fr with the following reference IE39108

All characteristics are given for information only and may be changed