

BAKES DUCT DILATOR 11 MM

Réf. IE39111



BAKES DUCT DILATOR 11 MM



Réf. IE39111

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040017631
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,2	Width packaging (cm):	
Length (cm) :	29,9	Length packaging (cm) :	
Height (cm) :	1,2	Height packaging (cm) :	
Weight (kg) :	0,048	Weight packaging (kg) :	0,048
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 21:06 - This article is available on our website : www.gsh-med.fr with the following reference IE39111

All characteristics are given for information only and may be changed