

DESJARDIN DUCT DILATOR N°5

Réf. IE40205



DESJARDIN DUCT DILATOR N°5



Réf. IE40205

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040017686
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,8	Width packaging (cm):	
Length (cm) :	30,0	Length packaging (cm) :	
Height (cm) :	0,8	Height packaging (cm) :	
Weight (kg) :	0,034	Weight packaging (kg) :	0,034
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 21:06 - This article is available on our website : www.gsh-med.fr with the following reference IE40205

All characteristics are given for information only and may be changed