

VOLKMANN RETRACTOR-SHARP 5 PRONGS

Réf. IE43521



VOLKMANN RETRACTOR-SHARP 5 PRONGS



Réf. IE43521

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040017839
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,3	Width packaging (cm):	
Length (cm) :	22,0	Length packaging (cm) :	
Height (cm) :	1,4	Height packaging (cm) :	
Weight (kg) :	0,098	Weight packaging (kg) :	0,098
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:12 - This article is available on our website : www.gsh-med.fr with the following reference IE43521

All characteristics are given for information only and may be changed