

RETRACTOR AUVERT (ABDOMINAL) - 260 X 340 MM

Réf. IE90401



RETRACTOR AUVERT (ABDOMINAL) - 260 X 340 MM



Réf. IE90401

FEATURES			
Medical device	No		
Marquage	Barcode	3661040178066	
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	1,105	Weight packaging (kg) :	1,105
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:35 - This article is available on our website : www.gsh-med.fr with the following reference IE90401

All characteristics are given for information only and may be changed