

DIRECTOR AND TONGUE TIE 16 CM

Réf. IS02416



DIRECTOR AND TONGUE TIE 16 CM



Réf. IS02416

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030050
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,0	Width packaging (cm):	
Length (cm) :	16,0	Length packaging (cm) :	
Height (cm) :	0,6	Height packaging (cm) :	
Weight (kg) :	0,009	Weight packaging (kg) :	0,009
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:38 - This article is available on our website : www.gsh-med.fr with the following reference IS02416

All characteristics are given for information only and may be changed