


NOVACK UTERINE SOUND, 23 CM

Réf. IS03300



NOVACK UTERINE SOUND, 23 CM



Réf. IS03300

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030104
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,7	Width packaging (cm):	
Length (cm) :	23,1	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,019	Weight packaging (kg) :	0,019
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:38 - This article is available on our website : www.gsh-med.fr with the following reference IS03300

All characteristics are given for information only and may be changed