

METAL FEMALE CATHETER 15 CM

Réf. IS03402



METAL FEMALE CATHETER 15 CM



Réf. IS03402

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030128
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,1	Width packaging (cm):	
Length (cm) :	16,2	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,012	Weight packaging (kg) :	0,012
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:39 - This article is available on our website : www.gsh-med.fr with the following reference IS03402

All characteristics are given for information only and may be changed