

METAL FEMALE CATHETER 15 CM (SET OF 4 PCS DIAM 2,3,4,5 MM)

Réf. IS03405



METAL FEMALE CATHETER 15 CM (SET OF 4 PCS DIAM 2,3,4,5 MM)



Réf. IS03405

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe Ir | |
| Marquage | CE | Barcode | 3661040087597 |
| Guarantee (month): | 24 | Customs Code | 9018908500 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 1,3 | Width packaging (cm): | |
| Length (cm) : | 16,3 | Length packaging (cm) : | |
| Height (cm) : | 0,7 | Height packaging (cm) : | |
| Weight (kg) : | 0,011 | Weight packaging (kg) : | 0,011 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 21:39 - This article is available on our website : www.gsh-med.fr with the following reference IS03405

All characteristics are given for information only and may be changed