

COLLIN VAGINAL SPECULUM 35 MM X 110 MM

Réf. IS04635



COLLIN VAGINAL SPECULUM 35 MM X 110 MM



Réf. IS04635

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030234
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,5	Width packaging (cm):	
Length (cm) :	11,0	Length packaging (cm) :	
Height (cm) :	6,3	Height packaging (cm) :	
Weight (kg) :	0,256	Weight packaging (kg) :	0,256
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:39 - This article is available on our website : www.gsh-med.fr with the following reference IS04635

All characteristics are given for information only and may be changed