

COLLIN VAGINAL SPECULUM 38 MM X 110 MM

Réf. IS04638



COLLIN VAGINAL SPECULUM 38 MM X 110 MM



Réf. IS04638

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040030241
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,8	Width packaging (cm):	
Length (cm) :	11,0	Length packaging (cm) :	
Height (cm) :	6,0	Height packaging (cm) :	
Weight (kg) :	0,274	Weight packaging (kg) :	0,274
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:39 - This article is available on our website : www.gsh-med.fr with the following reference IS04638

All characteristics are given for information only and may be changed