

## COLLIN VAGINAL SPECULUM 40 MM X 110 MM

Réf. IS04640



# COLLIN VAGINAL SPECULUM 40 MM X 110 MM



Réf. IS04640

FEATURES			
Medical device	Yes	Classe I r	
Marquage	CE	Barcode	3661040030258
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	4,0	Width packaging (cm):	
Length (cm) :	11,0	Length packaging (cm) :	
Height (cm) :	6,7	Height packaging (cm) :	
Weight (kg) :	0,280	Weight packaging (kg) :	0,290
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:39 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IS04640

All characteristics are given for information only and may be changed